

## SPINERGY SPOX 25" IS REPLACED BY SPINERGY X 25"



SPINERGY SPOX



SPINERGY X

### TRANSITION FROM REAR WHEEL SPINERGY SPOX TO SPINERGY X.

Shortly we will start to replace rear wheel Spinerger Spox 25" with Spinerger X 25".

Spinerger X is the standard rear wheel for Panthera X and it's a wheel with very low weight, high quality and stability.

The attachment of the push rim has the same design as all our other rear wheels which facilitate store-keeping. All our push rims fits all our rear wheels. Detachment and mounting of the pushrim is made from the outside, i.e. you don't need to remove the tire and tube for this operation.

All sizes of the Spox-rear wheel will eventually be replaced by Spinerger X rear wheels.

The price for Spinerger X will be exactly the same as for all versions of Spinerger Spox.

Example of article numbers

art.no. 1100130 Rearwheel 25" SPOX std, without pushrim, Schwalbe Right Run tire

replaced by

art.nor. 1100155 Rearwheel 25" Spinerger X, without pushrim, Schwalbe Right Run tire