

# panthera®

art.nr. mont-tippskydd-hojder Rev.2010-01

## Val av tippskydd vid olika förhöjningsatser

Förhöjning	Drivhjul	Tippskydd +0 mm	Tippskydd +30 mm	Tippskydd +60 mm
		art.nr. 4351000	art.nr. 4351030	art.nr. 4351060
+ 0 mm	24"	X		
+ 0 mm	25"	X		
+ 0 mm	26"	X		
+ 10 mm	24"	X		
+ 10 mm	25"	X		
+ 10 mm	26"	X		
+ 20 mm	24"	X		
+ 20 mm	25"	X		
+ 20 mm	26"	X		
+ 30 mm	24"	X		
+ 30 mm	25"		X	
+ 30 mm	26"		X	
+ 40 mm	24"		X	
+ 40 mm	25"		X	
+ 40 mm	26"		X	
+ 50 mm	24"		X	
+ 50 mm	25"		X	
+ 50 mm	26"		X	
+ 60 mm	24"			X
+ 60 mm	25"			X
+ 60 mm	26"			X
+ 70 mm	24"			X
+ 70 mm	25"			X
+ 70 mm	26"			X

# panthera®

## Choice of anti-tip device at various raising-kit

Raising-kit	Rearwheel	Anti-tip +0 mm	Anti-tip +30 mm	Anti-tip +60 mm
		art.nr. 4351000	art.nr. 4351030	art.nr. 4351060
+ 0 mm	24"	X		
+ 0 mm	25"	X		
+ 0 mm	26"	X		
+ 10 mm	24"	X		
+ 10 mm	25"	X		
+ 10 mm	26"	X		
+ 20 mm	24"	X		
+ 20 mm	25"	X		
+ 20 mm	26"	X		
+ 30 mm	24"	X		
+ 30 mm	25"		X	
+ 30 mm	26"		X	
+ 40 mm	24"		X	
+ 40 mm	25"		X	
+ 40 mm	26"		X	
+ 50 mm	24"		X	
+ 50 mm	25"		X	
+ 50 mm	26"		X	
+ 60 mm	24"			X
+ 60 mm	25"			X
+ 60 mm	26"			X
+ 70 mm	24"			X
+ 70 mm	25"			X
+ 70 mm	26"			X