

Adjusted by:

E-mail:

Adjusted for:

Tel:

Date:

Total weight (g):

The weight changes based on the choices made below.

S3 short low has a adjustable rear axle in carbon fibre. To find a balance that suits you the rear axle can be adjusted into different positions. Note desired measure at "Settings". S3 short has a seat height of 47 cm in the front and 43 at the rear, without cushion. To adjust the seat height you can chose other rearwheel sizes or a raising Kit with suitable castor fork from the List of options. **S3 short low** is 50 mm shorter and 25 mm lower than the S3.

Max userweight: 100 kg for seat width 33-39

Save the guide to the "desktop" and open it with [Adobe Reader](#), for best function Design your wheelchair by clicking in the circles in this guide. Blue, underlined text is clickable to [show info](#) or image. BB = seatwidth, HH= backrest height

Original design

Alternative options

Other options

[info](#)

1) Chassis-colour

- Black** - Anodic Black
- Clear

Optional colour - A higher price and longer delivery time for optional colour, Fili in art.no. in the square

2) Seat width (BB) [article number](#)

- 33 cm
- 36 cm
- 39 cm

Art.no include: chassis, seat, calf strap, backrest lock

3) Seat depth

- Clear
- Seat depth 35 cm**
Std for all seat width

- Adjustable seat depth 27-33 cm**
For seat width 33-39 cm

4) Backrest

- Standard**
523BBHH VER2
- Clear

Height (HH)

- 20 cm
- 25 cm
- 30 cm
- 35 cm
- 40 cm
- 45 cm
- Clear

Tapered 525BBHH VER2
Curved back 524BBHH VER2
Info With backward curving the anti-tip device is complemented with "Knob anti-tip device S3 low profile", art.no. 4353121

Other backrest

art. _____

5) Rear axle

- Carbon standard**
33023BB
- Clear

Carbon +15 mm
33023BC* * seat width BB +1
Use this rear axle when choosing arm support at Accessories and/or increased seat height

Other rear axle or seat height

art. _____

6) Castor and fork

- 120 mm [art.no.](#)
- Clear

- 90 mm [art.no.](#)
- 150 mm [art.no.](#)

7) Footrest

- Footrest Titanium**
36001BB
- Footrest U-shape small heel bar**
3600134 (width 33)
- Clear

- Titanium, extended, reinforced**
36111 BB needed for Freewheel
- Foldable footplates**
[art.no.](#)
- Heelstrap for footplates**
[art.no.](#)
- Plastic plate**
36010BB

Other footrest

art. _____

8) Brake

- Brake S3**
4630000
- Clear

- One hand brake right**
[art.no.](#)

Other brake

art. _____

9) Rear wheel - without pushrim

- 24" **standard**
[art.no.](#)
- Clear

- 22" **standard** Rear wheel sizes
[art.no.](#)
- 25" **Spinergy X** black spokes
[art.no.](#)
- 26" **Spinergy X** black spokes
[art.no.](#)

Other rear wheel

art. _____

Tire - choose tire

- Schwalbe Right Run** 22",24",25",26"
- Clear

- Rugged** 24"
- Schwalbe **Marathon Plus** 24",25",26"

PU 24" (massive)

reTyre-tyre skins (does not fit PU tyres)

10) Pushrim

- 24" Titanium**
190DD11 (x2 / wheelchair)
- Clear DD = Wheel diameter

- Aluminum**
1904404 x2 (04=24", 05=25", 06=26")
- Paragrip**
19051DDBT x2
- Tetragrip**
19052DDBT x2

Maxgrip 22", 24" (contains latex)
1903060 x2 (60=24", 80=22")



Accessories



Side guard

- [Side guard S3](#) [Video](#)
[art.no.](#)

With 24" rearwheels,
choose Side guard 26"

- Clear



Arm support

- [Arm support S3](#) [Video](#)
[art.no.](#)

With 24" rearwheels,
choose Arm support 26"



Side guard Carbon

- [Side guard Carbon](#)
[art.no.](#)



Push handles

- [Height adjustable S3](#) [Video](#)
4252060 VER2

- [Height adjustable S3 +10](#)
4252063 VER2

Max. userweight 100kg

- Clear



Push handle Foldable

- [Foldable S3](#) [Video](#)
4245000



Push bar

- [Push bar S3](#)
42400BB VER2

Max. userweight 100kg



Anti-tip low chassis

- [Anti-tip std seat height](#)
4355000 [Video](#)

[Other seat height](#)




Cushion Wedge cushion

- [2,5 cm](#)
572BB35

- [5 cm](#)
575BB35

- Clear

- [Wedge cushion](#)
57710BB



Spoke protectors

- [24" logo](#)
4500000

- [24" no logo](#)
4500011

- Clear

- [22" logo](#)
4500100

- [22" no logo](#)
4500111



Other

- [Tube protection front](#)

9900251 (2 pcs)

- [Tube protection rear](#)

9900251 (1 pcs)

- [Crutch holder S3](#)

9900301 (x1/crutch) _____ pcs

- [Backrest wedge](#)

4710012 _____ pcs

- [Table S3](#)

6810020



Chose art.no. for Freewheel according to the **footrest height**
of the wheelchair

Read more about Freewheel on the [homepage](#)

- [Freewheel](#)

Set-up adjustments



Back rest angle

Note measure (A) for desired back rest angle.
14-28 mm. Standard=20 mm

A= _____ mm



Rear axle adjustment

Note measure (A) for desired balance.
Standard=40 mm

A= _____ mm



Footrest height [Footrest - Foot plates](#)

Note measure (A) and (B) for desired footrest.
[Bambino foot plate with Feewheel - info](#)

A= _____ mm B= _____ mm



Transportation marking

- YES do mark the wheelchair with
 attachment points for transport in
disabled service vehicle.

Art.no.9000152

NOTE! See [manual](#) for more information.



Extra