

## Panthera S3 / Bambino 3

Framgaffel med länkhjul vid olika storlekar på drivhjul

<b>120 mm</b>			
Drivhjul	Framgaffel	art.nr.	
22"		<b>2100020</b> Gaffel 75, hjul 120, tapp 3	x2
24"		<b>2100070</b> Gaffel 105, hjul 120, tapp 3	x2
25"		<b>2100100</b> Gaffel 120, hjul 120, tapp 3	x2
26"		<b>2100120</b> Gaffel 105, tapp 33	x2
		<b>2002020</b> Länkhjul S3, 120	x2
<b>90 mm</b>			
Drivhjul	Framgaffel	art.nr.	
22"		<b>2100030</b> Gaffel 90, hjul 90, tapp 3	x2
24"		<b>2100090</b> Gaffel 120, hjul 90, tapp 3	x2
25"		<b>2100120</b> Gaffel 105, tapp 33	x2
		<b>2002010</b> Länkhjul S3, 90	x2
26"		<b>2100130</b> Gaffel 120, tapp 33	x2
		<b>2002010</b> Länkhjul S3, 90	x2
<b>150 mm</b>			
Drivhjul	Framgaffel	art.nr.	
22"			
24"		<b>2100050</b> Gaffel 90, hjul 150, tapp 3	x2
25"		<b>2100080</b> Gaffel 105, hjul 150, tapp 3	x2
26"		<b>2100110</b> Gaffel 120, hjul 150, tapp 3	x2

## ENG

Castor fork with castor for different rearwheel sizes.

<b>120 mm</b>			
Rearwheel	castro fork	art.no.	
22"		<b>2100020</b> Fork 75 Castor 120 Bolt 3	x2
24"		<b>2100070</b> Fork 105 Castor 120 Bolt 3	x2
25"		<b>2100100</b> Fork 120 Castor 120 Bolt 3	x2
26"		<b>2100120</b> Fork 105 Bolt 33	x2
		<b>2002020</b> Frontcaster S3 120 mm	x2
<b>90 mm</b>			
Rearwheel	castro fork	art.no.	
22"		<b>2100030</b> Fork 90 Castor 90 Bolt 3	x2
24"		<b>2100090</b> Fork 120 Castor 90 Bolt 3	x2
25"		<b>2100120</b> Gaffel 105, tapp 33	x2
		<b>2002010</b> Frontcaster S3 90 mm	x2

26"	<b>2100130</b>	Fork 120 Bolt 33	x2
	<b>2002010</b>	Frontcaster S3 90 mm	x2
	<b>150 mm</b>		

Rearwheel castro fork  
**art.no.**

22"			
24"	<b>2100050</b>	Fork 90 Castor 150 Bolt 3	x2
25"	<b>2100080</b>	Fork 105 Castor 150 Bolt 3	x2