


Eploded View Panthera

Footplate Carbon Angle-adjustable

permobil



Pos.	Benämning		Art.nr	Antal
1	Plate carbon	Not sold separately		1
2	Attachment foot plate carbon left	Not sold separately	6506011	1
3	Attachment foot plate carbon right	Not sold separately	6506012	1
4a	Bolt BN 1206 M5x12 Oxidized Black	Not sold separately	6004512C	4
4b	Bolt BN 1206 M5x16 Oxidized Black	Not sold separately	6004516C	4
5	Bolt MELC6S M5x8 FZSV	Not sold separately	6004588	4
6	Washer M5x13x1	Not sold separately	8000513	4
7	Lock nut DIN 985 steel M5 fzb		7000050	4
8	Attachment angle-adjust. foot plate		6506050	2
9	Clamp B3 Side cpl. right		4234042	1
10	Clamp B3 Side cpl. left		4234041	1
11	Bolt MC6S M6x18 12.9 FZSV		6002648	4
12	Calf strap Panthera			1