


X BACK REST UPHOLSTERY UPPER PART

Click on images for larger format



UPPER PART BACK REST UPHOLSTERY

THE UPPER PART FOR X BACK REST UPHOLSTERY IS REPLACED WITH S3 UPHOLSTERY UPPER PART

Below you can see the article numbers that replaces the upper part on the X back rest upholstery. For back rest height 20-24 cm, use article numbers in the left column and for height 25-36 cm use article numbers in the right column. The article numbers is complete with both left and right side.

		<i>One article number for both left & right side</i>			
EXPIRE		REPLACES (back rest height 20-24 cm)		REPLACES (back rest height 25-36 cm)	
ART.NO.	DESCRIPTION	ART.NO.	DESCRIPTION	ART.NO.	DESCRIPTION
3266433	Backupholstery X upper left w 33	5266233	Upholst. backrest upperpart w 33-20, with strap	5266133	Upholstery backrest upperpart w 33
3266533	Backupholstery X upper right w 33				
3266436	Backupholstery X upper left w 36	5266236	Upholst. backrest upperpart w 36-20, with strap	5266136	Upholstery backrest upperpart w 36
3266536	Backupholstery X upper right w 36				
3266439	Backupholstery X upper left w 39	5266239	Upholst. backrest upperpart w 39-20, with strap	5266139	Upholstery backrest upperpart w 39
3266539	Backupholstery X upper right w 39				
3266442	Backupholstery X upper left w 42	5266242	Upholst. backrest upperpart w 42-20, with strap	5266142	Upholstery backrest upperpart w 42
3266542	Backupholstery X upper right w 42				
3266445	Backupholstery X upper left w 45	5266245	Upholst. backrest upperpart w 45-20, with strap	5266145	Upholstery backrest upperpart w 45
3266545	Backupholstery X upper right w 45				