

## Product news Push rims with PRC and GP on SPOX

All SPOX dimensions will now be available in PRC (grey vinyl) and GP (gel-coated plastic). They are attached to the wheel rims at 6 points. NB – equal distance of 38 mm (centre wheel rim to centre push rim) between the wheel rims and push rims. Non-adjustable. This provides an air margin of 12 mm on PRC and 9 mm on GP (between the wheel rim and push rim).



Spox with PRC



Spox with GP

### Product

Rear wheel 24" Spox standard PRC rear wheel  
Rear wheel 24" Spox standard gel plastic rear wheel  
Rear wheel 24" Spox standard PRC rear wheel red  
Rear wheel 24" Spox standard gel plastic rear wheel red  
Rear wheel 25" Spox standard PRC rear wheel  
Rear wheel 25" Spox standard gel plastic rear wheel  
Rear wheel 26" Spox standard PRC rear wheel  
Rear wheel 26" Spox standard gel plastic rear wheel

### Article number

110 00 85  
110 00 86  
110 00 87  
110 00 88  
110 01 33  
110 01 32  
110 01 04  
110 01 05