

# panthera®

Info sheet 2011-01

2011-04-04

## Quick release for onearm drive

From 1st of January 2011 we deliver all onearm drive with quick release (QR).

This makes it easier to remove the rearwheels when transporting the wheelchair.

After removing the spring loaded drive shaft (3) and loosen two locking pins and you can easily pull off the rearwheels from their attachments.

Onearm drive is used together with the straight, steel rearaxle (2) without camber, that is ordered separately as for the drive shaft to get a complete onearm drive kit.

When ordering onearm drive please choose size of the rearwheels and what kind of pushrim, aluminum or friction. Also choose steel rear axle and drive shaft for the current seat width of the chair.

Onearm drive fits all wheelchairs except Micro, BT, U2 light and Panthera X.

See next page for article numbers!



# panthera®



## Article numbers for onearm drive with QR:

Contains 2 rearwheels with attachment and 2 pushrims mounted on one rearwheel.

art.no. 1101520	Rearwheel 20" onearm drive with QR	with aluminum-rim pair complete	1
art.no. 1101530	Rearwheel 20" onearm drive with QR	with friction-rim pair complete	
art.no. 1101522	Rearwheel 22" onearm drive with QR	with aluminum-rim pair complete	1
art.no. 1101532	Rearwheel 22" onearm drive with QR	with friction-rim pair complete	
art.no. 1101524	Rearwheel 24" onearm drive with QR	with aluminum-rim pair complete	1
art.no. 1101534	Rearwheel 24" onearm drive with QR	with friction-rim pair complete	
art.no. 1101526	Rearwheel 26" onearm drive with QR	with aluminum-rim pair complete	1
art.no. 1101536	Rearwheel 24" onearm drive with QR	with friction-rim pair complete	

## Article number for steel rearaxle:

art.no. 3301127	Tube rearaxle w bult 27	(chair width 27)	2
art.no. 3301130	Tube rearaxle w bult 30	(chair width 30)	
art.no. 3301133	Tube rearaxle w bult 33	(chair width 33)	2
art.no. 3301136	Tube rearaxle w bult 36	(chair width 36)	
art.no. 3301139	Tube rearaxle w bult 39	(chair width 39)	2
art.no. 3301142	Tube rearaxle w bult 42	(chair width 42)	
art.no. 3301145	Tube rearaxle w bult 45	(chair width 45)	2
art.no. 3301150	Tube rearaxle w bult 50	(chair width 50)	

## Article number for drive shaft:

art.no. 1100360	Drive-shaft one-armdrive conic 195-235/27	(chair width 27)	3
art.no. 1100361	Drive-shaft one-armdrive conic 220-275/30-33	(chair width 30-33)	
art.no. 1100362	Drive-shaft one-armdrive conic 255-315/33-36	(chair width 33-36)	3
art.no. 1100363	Drive-shaft one-armdrive conic 295-370/39-42	(chair width 39-42)	
art.no. 1100364	Drive-shaft one-armdrive conic 350-460/42-50	(chair width 42-45)	3