

panthera®

Info sheet 2011-05

2011-09-20

Changes on U2 light

We will gradually introduce changes to the Panthera U2 light during a transitional period. The changes are a new design and some technical changes. However, with the same price as before.

Changes are; new stickers with the Panthera logotype, new design of the calf band and back rest upholstery. The caster fork has a new design. The back rest attachment has a new design.

- 1) New stickers with the Panthera logotype and new design of the calf band gives a nice and modern look.
- 2) The castor fork has a new design.
- 3) The back rest attachment has a new design and the attachment for the side-guard is lowered to fit the new Carbon sideguards..
- 4) The back rest upholstery is made of the same material as you find on the Panthera X and due to it's structure, airy and comfortably to sit against.

New art.no:

Fork U2 light X cpl w. std 3" wheel art.no. 2300180

Backrest cpl width 33 height 25 light B art.no. 3233325LB
+ all backrest for chair width between 33 cm and 45 cm and heights between 25 cm and 35 cm

Backrest upholstery width 33 height 25 light B art.no. 3263325LB
+ all backrest upholstery for chair width between 33 cm and 45 cm and heights between 25 cm and 35 cm



This information is available at our webb: [www.panthera.se / english / info / technical info](http://www.panthera.se/english/info/technical%20info)

Panthera AB, Gunnebogatan 26, SE-163 53 Spånga, Sweden, +46-8-761 50 40, www.panthera.se