

panthera®

Info sheet 2011-06

2011-09-12

Panthera X brake



Now we have the new Panthera X brake ready for delivery. Chairs delivered as of the 12th September 2011 will have the new brake mounted, if ordered.

The new brake weigh 258 gram, which is 150g less than the old brake.

The brake is made of aluminum and has parts of carbon fiber.

It has improved geometry and accuracy in the brake movement which make it is easier to use for those who have reduced function in the hand.

The design and location of the linkrod is not in the way when lifting the chair.



Description

art.no.

Brake onehand X 33 Right	4623302
Brake onehand X 33 Left	4623301
Brake onehand X 36 Right	4623602
Brake onehand X 36 Left	4623601
Brake onehand X 39 Right	4623902
Brake onehand X 39 Left	4623901
Brake onehand X 42 Right	4624202

PANTHERA TECHNOLOGY

This info is also available at our webb: [www.panthera.se / english / info / technical-info](http://www.panthera.se/english/info/technical-info).

Panthera AB, Gunnebogatan 26, SE-163 53 Spånga, Sweden, +46-8-761 50 40, www.panthera.se