

panthera®

Info sheet 2011-08

2011-10-24

Bambino with abduction front.



Now it is possible to order Bambino with abduction front, which means that the frame tubes by the knees are slightly flared to give more room. The width between the frame tubes has increased with 6 cm. Also the footrest is 6 cm wider. When ordering a footrest for an abduction front then it should be 6 cm wider than the seatwidth. Width 33 cm has 6 weeks delivery time.

New article numbers:

Art.no. 3572400G	width 24 cm
Art.no. 3572700G	width 27 cm
Art.no. 3573000G	width 30 cm
Art.no. 3573300G	width 33 cm

