

# panthera®

Info sheet 2011-09

2012-03-19

## Freewheel



Panthera can now offer its customers a new accessory, an invention called Freewheel. It is developed in the United States by an American named Pat Dougherty.

Freewheel - is an accessory that is easily mounted on the wheelchair's footrest as shown in the image. It dramatically increases accessibility for the wheelchair on surfaces such as lawns, snow-covered parking lots, paths, gardens, curbs, etc.

It is intended for use outdoors. The wheelchair's length increases by 46 cm when Freewheel is mounted. The assembly and disassembly of Freewheel is very easy and can be done by the user himself sitting in the wheelchair.

An attachment to hang the Freewheel on the backrest of the wheelchair when it's not in use, is also available.

See video on Freewheel on Panthera's homepage. [Link](#)

**Artikelnummer:**

**Art.no. 2600102**

