


panthera®

Info sheet 2012-01

2013-03-25

New standard hub



Panthera have a new supplier of hubs.

We want this action to ensure a high level of quality and reliability.

Panthera will gradually introduce the new hub on our standard rearwheels, starting from January 1, 2012.

The new hub features a new design and a different alloy of aluminum. Also the axle has a new design. The new rearwheel axle is backward compatible with the old hub and the old rearwheel axle will fit the new hub.

The new hub is easily recognized by the black color and wheel axle of the beveled button and the laser-written Panthera logo.

Hubshell with bearings art.no. 1700070
Axle QR 1/2" / 12,7mm x 91 mm SJ art.no. 1100022