

panthera®

Info sheet 2012-03

2012-02-13

Assistant brake



Panthera can now offer an assistant brake kit .

The kit contains two 24" rear wheels with drum brakes and aluminum push rims and two brake handles with wires. Also containing assembly parts for mounting the rear wheels on the old std rear axle and for mounting the brake handles on the height adjustable push handles / pushbar.

If the chair is not equipped with the old std rear axle and height adjustable push handles / pushbar the order needs to be completed with those articles.

The rear wheels are 24" as standard but can be ordered as 20" or 22" but they are not stock-goods and has longer delivery time.

Assistant brake fits:

Panthera S2, S2 short, S2 short-low, S2 Large, Swing, Swing-short, U2 and Bambino

Rear wheel 24" with assistant break incl Alu push rims

art.nr. 1103241

Rear wheel 24" with assistant break without pushrim

art.nr. 1103240