

panthera®

Info sheet 2012-04

2012-05-11

Widths of the Bambino foot plate



The widths of the new Bambino foot plate.

The new Bambino foot plate is 24 mm wider than before. This is to better use the width of the chassis and create more space for the feet. Previously, you could order a foot plate, one size larger, to get more space for the feet, but that is not possible now. You use the same width on the foot plate as the width of the chassis.

However with **Bambino with abduction front** you use these widths:

Panthera Bambino 24 w. abduction front	Footplate Bambino 30 cm +24 mm
Panthera Bambino 27 w. abduction front	Footplate Bambino 33 cm +24 mm
Panthera Bambino 30 w. abduction front	Footplate Bambino 36 cm +24 mm
Panthera Bambino 33 w. abduction front	Footplate Bambino 39 cm +24 mm

When using the **Bambino foot plate on Panthera U2** you choose a smaller size than the chassis width, see table:

Panthera U2 33 cm	Footplate Bambino 24 cm +24 mm
Panthera U2 36 cm	Footplate Bambino 27 cm +24 mm
Panthera U2 39 cm	Footplate Bambino 30 cm +24 mm
Panthera U2 42 cm	Footplate Bambino 33 cm +24 mm
Panthera U2 45 cm	Footplate Bambino 36 cm +24 mm

The new Bambino footplates:

Footplate Bambino 24 cpl kit+24mm	art.no. 4002125
Footplate Bambino 27 cpl kit+24mm	art.no. 4002128
Footplate Bambino 30 cpl kit+24mm	art.no. 4002131
Footplate Bambino 33 cpl kit+24mm	art.no. 4002134
Footplate Bambino 36 cpl kit+24mm	art.no. 4002137
Footplate Bambino 39 cpl kit+24mm	art.nr. 4002140