



Info sheet 2013-01

15/03/2013

New 24" standard rear wheel

Beginning in March 2013, will gradually begin to deliver our wheelchairs with a new and improved, 24" standard rear wheel. The improvements entail an even stronger and more rigid rim, as well as a new attachment for the push rim with 6 attachment points that provide the capability to change the push rim without needing to remove the tyre and tube, since the screws are accessible from the outside of the rim.

Moreover, all push rim types will fit one and the same wheel; titanium, aluminium PRC, gel coated and maxgrip push rims will all have the same attachment and fit the new wheel.

The hub is the same as previously updated in January 2012. This entails that the quick-release axle is the same as before. The new wheels fit all chairs that the old wheel fits on.

Older push rims that fit the older standard wheel will not fit the new since they have different inside diameters on the rim. However, push rims that fit the Spinergy wheels will fit the new standard rear wheel.

The new wheels fit the S2, S2 short, S2 short-low, S2 large, S2 swing, S2 swing-short, U2 and Bambino.

New article numbers:

art.no. 1300000	Rear wheel 24" 6P without push rim
art.no. 1300000R	Rear wheel 24" 6P without push rim without tyre
art.no. 1300001	Rear wheel 24" 6P without push rim PU tyre
art.no. 1300010	Rear wheel 24" 6P Push rim Alu
art.no. 1300010R	Rear wheel 24" 6P Push rim Alu without tyre
art.no. 1300011	Rear wheel 24" 6P Push rim Alu PU tyre
art.no. 1300020	Rear wheel 24" 6P Push rim Titan
art.no. 1300020R	Rear wheel 24" 6P Push rim Titan without tyre
art.no. 1300021	Rear wheel 24" 6P Push rim Titan PU tyre
art.no. 1300030	Rear wheel 24" 6P Push rim MaxGr thin
art.no. 1300030R	Rear wheel 24" 6P Push rim MaxGr thin without tyre
art.no. 1300031	Rear wheel 24" 6P Push rim MG PU tyre
art.no. 1300040	Rear wheel 24" 6P Push rim Gel coated
art.no. 1300040R	Rear wheel 24" 6P Push rim Gel coated without tyre
art.no. 1300041	Rear wheel 24" 6P Push rim GP PU tyre
art.no. 1300050	Rear wheel 24" 6P Push rim PRC gray
art.no. 1300050R	Rear wheel 24" 6P Push rim PRC gray without tyre
art.no. 1300051	Rear wheel 24" 6P Push rim PRC PU tyre

