



Info sheet 2013-03
2013-02-28

Push rims PRC and foam coated expires

As from the 1st of March 2013 we will **gradually** replace the push rims PRC and foam coated with new models. This goes for all sizes.

The reason for replacing these push rims is that our American supplier has unexpectedly closed their factory and stopped all deliveries.

Therefore we need to introduce the new push rims, that Panthera has tested and evaluated during the autumn / winter 2012 / 2013, with a good result, earlier than expected.

It is the push rims Tetra Grip and Para Grip, from the well-reputed German manufacturer, Tetra Equipment.

PRC (gray color) replaced by **Para Grip** (black color)*
Foam coated (soft) replaced by **Tetra Grip** (slightly harder)*

*These push rims are mounted with a M5 screw instead of the screw with a thread in inch.

PRC	art.no.	replaced by	Para Grip	art.no.
Push rim PRC for 24"-wheel	1903300	" "	Push rim Para 24" 4P C/C 507	1903406
Push rim PRC 16" short dist spokes	1903303	" "	Push rim Para 20" 6 tabs	1905120
Push rim PRC 18" short dist spokes	1903305	" "	Push rim Para 22" 6 tabs	1905122
Push rim PRC for 24" SPOX 6P	1904304	" "	Push rim Para 24" 6 tabs	1905124
Push rim PRC for 25" SPOX	1904305	" "	Push rim Para 25" 6 tabs	1905125
Push rim PRC for 26" SPOX	1904306	" "	Push rim Para 26" 6 tabs	1905126
Foam coated			Tetra Grip	
Push rim Foam coated for 24"-wheel	1903200	" "	Push rim Tetra 24" 4P C/C 507	1905324
Push rim Foam coated for 24" SPOX	1904204	" "	Push rim Tetra 24" 6 tabs	1905224
Push rim Foam coated for 25" SPOX	1904205	" "	Push rim Tetra 25" 6 tabs	1905225
Push rim Foam coated for 26" SPOX	1904206	" "	Push rim Tetra 26" 6 tabs	1905226