

New standard rear axle for Bambino



Soon, the Bambino will be delivered with a new standard rear axle.

The new rear axle has a camber angle of 6° and the tube, made of aluminum, has a diameter of 26 mm.

The advantages of the new rear axle model is improved driving characteristics and that it fits all variations of the Bambino chassis.

art.no.	Description	Fits
Complete kit		
3362416	Camber axle kit 24 6 dgr 26 mm	Bambino std + Bambino Abduction width 24 cm
3362716	Camber axle kit 27 6 dgr 26 mm	Bambino std + Bambino Abduction width 27 cm
3363016	Camber axle kit 30 6 dgr 26 mm	Bambino std + Bambino Abduction width 30 cm
3363316	Camber axle kit 33 6 dgr 26 mm	Bambino std + Bambino Abduction width 33 cm
Rear axle tube		
3362406	Camber axle tube 26 mm/24 cm 6 gr kpl	Bambino std + Bambino Abduction width 24 cm
3362706	Camber axle tube 26 mm/27 cm 6 gr kpl	Bambino std + Bambino Abduction width 27 cm
3363006	Camber axle tube 26 mm/30 cm 6 gr kpl	Bambino std + Bambino Abduction width 30 cm
3363306	Camber axle tube 26 mm/33 cm 6 gr kpl	Bambino std + Bambino Abduction width 33 cm
Attachment		
3302510	Assembly kit Carbon axle +10mm	Bambino std + Bambino Abduction all widths