

## Overview of new rear wheels and push rims



New standard rear wheels.



Para Grip push rim supersedes PRC.



Tetra Grip push rim supersedes Foam coated.

### New standard rear wheels

As previously announced, our existing standard rear wheels are being superseded by new, improved ones featuring stronger, stiffer rims and six new attachment points for the push rim. This makes it possible to replace push rim without removing the tyre and tube since the screws are accessible from the outer side of the rim. The introduction began in March 2013.

### New replacements for PRC and foam coated push rims

Early in 2013, we were informed that our American supplier had stopped producing and delivering PRC and foam coated push rims, which took effect immediately. This forced us to bring forward the introduction of the new push rims, Tetra Grip and Para Grip from the well-known German manufacturer Tetra Equipment, which Panthera had evaluated and approved during autumn 2012/winter 2013. The Tetra Grip is somewhat harder than the the foam coated push rim and offers an excellent non-slip surface and greater durability. In terms of function, the Para Grip is very similar to the PRC but it is black.

### The following is a comprehensive list of our rear wheels and push rims

#### New standard rear wheels

art.no.	Description	supersedes
1320000	Rear wheel 20" 6P without push rim	1102001
1320020	Rear wheel 20" 6P push rim titanium	1102000
1320030	Rear wheel 20" 6P push rim MaxGrip thin	1100008
1310000	Rear wheel 22" 6P without push rim	1102201
1310020	Rear wheel 22" 6P Push rim titanium	1102200
1310021	Rear wheel 22" 6P push rim titanium PU tyre	1102202
1310030	Rear wheel 22" 6P push rim MaxGrip	1100009
1300000	Rear wheel 24" 6P without push rim	1100060
1300010	Rear wheel 24" 6P push rim aluminium	1100039
1300011	Rear wheel 24" 6P push rim aluminium PU tyre	
1300020	Rear wheel 24" 6P push rim titanium	1000000
1300021	Rear wheel 24" 6P push rim titanium PU tyre	1100005
1300030	Rear wheel 24" 6P push rim MaxGrip thin	1100001
1300031	Rear wheel 24" 6P push rim MaxGrip PU tyre	1102472
1300040	Rear wheel 24" 6P push rim Tetragrip	1100002
1300041	Rear wheel 24" 6P push rim Tetragrip PU tyre	1102471
1300050	Rear wheel 24" 6P push rim Paragrip	1100003
1300051	Rear wheel 24" 6P push rim Paragrip PU tyre	1102470

#### New push rims for the new standard rear wheels

art no.	Description
1905220	Push rim Tetragrip 20" 6P
1905120	Push rim Paragrip 20" 6P
1902011	Push rim Titanium 20" 6P 421 mm
1905222	Push rim Tetragrip 22" 6P
1905122	Push rim Paragrip 22" 6P
1902211	Push rim Titanium 22" 6P 459 mm
1905224	Push rim Tetragrip 24" 6P
1905124	Push rim Paragrip 24" 6P
1902411	Push rim Titanium 24" 6P 510 mm
1903060	Push rim MaxGrip 24" wheel 6P 510 mm

See the following page  
for more information:

Push rims		Superseded by:		Fits wheel	Fits wheel	Fits wheel	Fits wheel	Fits wheel	Rim-mounted	Spoke-mounted
art no.	Description	art no.	Description	20"	22"	24"	25"	26"		
1902400	Push rim titanium for 24" wheel	1902411	Push rim titanium for 24" wheel 6P 510 mm			x			x	
1902200	Push rim titanium for 22" wheel	1902211	Push rim titanium 22" 6P 459 mm		x				x	
1902000	Push rim titanium for 20" wheel	1902011	Push rim titanium 20" 6P 421 mm	x					x	
1903040	Push rim MaxGrip for 24" wheel	1903060	Push rim MaxGrip 24" wheel 6P 510 mm			x			x	
1903200	Push rim foam coate for 24" wheel	1905324	Push rim Tetra 24" 4P C/C 507			x*			x*	
1903300	Push rim PRC for 24" wheel	1903406	Push rim (Para) 24" 4P C/C 507			x*			x*	
1902511	Push rim titanium for 25" wheel	Same					x		x	
1904205	Push rim foam coated for 25" SPOX	1905225	Push rim Tetra 25" 6P				x		x	
1904305	Push rim PRC for 25" SPOX pc	1905125	Push rim Para 25" 6P				x		x	
1902611	Push rim titanium for 26" wheel Spox	Same						x	x	
1904206	Push rim foam coated for 26" SPOX pc	1905226	Push rim Tetra 26" 6P					x	x	
1904306	Push rim PRC for 26" SPOX pc	1905126	Push rim Para 26" 6P					x	x	
	New	1902430	Push rim ALU for 24" wheel 6P 510 mm			x			x	
1907020	Push rim MaxGrip 17" thin long spacer	Same		x	x	x	x	x		x
1907030	Push rim MaxGrip 17" thin short spacer	Same		x	x	x	x	x		x
1907320	Push rim MaxGrip 17" thick long spacer	Same		x	x	x	x	x		x
1907330	Push rim MaxGrip 17" thick short spacer	Same			x	x	x	x		x
1908020	Push rim MaxGrip 18" thin long spacer	Same			x	x	x	x		x
1908030	Push rim MaxGrip 18" thin short spacer	Same			x	x	x	x		x
1908320	Push rim MaxGrip 18" thick long spacer	Same			x	x	x	x		x
1908330	Push rim MaxGrip 18" thick short spacer	Same			x	x	x	x		x
1903303	Push rim PRC 16" short spacer pc spokes	Same		x**	x					x
1903305	Push rim PRC 18" short spacer pc spokes	1905122	Push rim Para 22" 6P						x	

\* For older 4-clip wheels \*\* Rim-mounted with special attachment