

FREEWHEEL - HIGH MOUNTED FOOTREST

The manufacturer now can offer variants of the Freewheel, adapted for high mounted footrest, where it before was impossible to use Freewheel. The measures state the height, in cm, between the floor and the top of the footrest.

The new variants, incl. std are:

art.no.	Freewheel	Lowest height	Highest height
2600102	Freewheel std		
2600105	Freewheel special height 5,5" (14 cm)	11,5 cm	16,5 cm
2600106	Freewheel special height 6,5" (16,5 cm)	14 cm	19 cm
2600107	Freewheel special height 7,5" (19 cm)	16,5 cm	21,5 cm
2600108	Freewheel special height 8,5" (21,5 cm)	19 cm	24 cm

Use the table below to determine which Freewheel to use together with your footrest.

NOTE! The combination of Freewheel and narrow footrest, width 30 or 33 cm, is not recommended because of the space for the feet then becomes too small when Freewheel is mounted. Contact Panthera for more information.



Model of wheelchair	Shims (A, B)	Dist. from the floor, in cm	Freewheel - model	Art.no.	Art.no.	Type of footrest
U2 / U3 / light / X	B	4	Std	2600102		Reinforced
U2 / U3 / light / X	B	6	Std	2600102		Reinforced
U2 / U3 / light / X	B	8	Std	2600102		Reinforced
U2 / U3 / light / X	B	10	Std	2600102		Reinforced
U2 / U3 / light / X	A (25" wheel)	11,2	Std	2600102		Reinforced
U2 / U3 / light / X	B (26" wheel)	12,5	5,5"	2600105		Reinforced
Model of wheelchair	Shims (A, B)	Dist. from the floor, in cm	Freewheel - model	Art.no.	Art.no.	Type of footrest
S2 / S3	B	4	Std	2600102		Reinforced
S2 / S3	B	6	Std	2600102		Reinforced
S2 / S3	B	8	Std	2600102		Reinforced
S2 / S3	B	10	Std	2600102		Reinforced
S2 / S3	A (25" wheel)	11,2	Std	2600102		Reinforced
S2 / S3	B (26" wheel)	12,5	5,5"	2600105		Reinforced
Model of wheelchair	Shims (A, B)	Dist. from the floor, in cm	Freewheel - model	Art.no.	Art.no.	Type of footrest
S2 / S3-short	B	4	-	-		-
S2 / S3-short	B	6	-	-		-
S2 / S3-short	B	8	Std with 3" castors	2600102		U-shape S2
S2 / S3-short	B	10	Std with 3" castors	2600102		U-shape S2
S2 / S3-short	A (25" wheel)	12	5,5"	2600105		U-shape S2
S2 / S3-short	B (26" wheel)	14	5,5"	2600105		U-shape S2
Model of wheelchair	Shims (A, B)	Dist. from the floor, in cm	Freewheel - model	Art.no.	Art.no. Bambino Kit *	Type of footrest
Bambino	Bambino kit*	8-14	5,5"	2600103	Art.no. 2600103 contains Bambino kit	Footplate
Bambino	Bambino kit	14-19	6,5"	2600106	+2602061	Footplate
Bambino	Bambino kit	16,5-21,5	7,5"	2600107	+2602061	Footplate
Bambino	Bambino kit	19-24	8,5"	2600108	+2602061	Footplate

* Bambino kit contains parts for adapting the footplate for Freewheel. (Footplate is not included in Bambino kit)