

PUSHRIM PARAGRIP AND TETRAGRIP WITH BENT TABS



NEW PARA- AND TETRAGRIP PUSHRIMS WITH BENT TABS



OLD PARA- AND TETRAGRIP PUSHRIMS WITH STRAIGHT TABS

Pushrim Paragrip and Tetragrip will soon be delivered with bent attachments (tabs), the same design as for the titanium pushrim. They will gradually replace the old with straight tabs.

The benefits will be larger outer diameter, 560 mm, previously 530 mm, which gives easier propelling of the wheelchair.

The new pushrims are mounted in the same way as the old and the measure between the tabs are the same.

To be able to separate the new pushrims from the old they are supplied with BT (bent tabs) in the article number.

You can still order the pushrims with straight tabs as long they are in stock.



DIFFERENCE IN DIAMETER BETWEEN THE OLD PUSHRIM AND THE NEW TO THE RIGHT - 530 TO 560 MM

art.no. 1905120BT	Pushrim Para 20" 6 bent tabs
art.no. 1905122BT	Pushrim Para 22" 6 bent tabs
art.no. 1905124BT	Pushrim Para 24" 6 bent tabs
art.no. 1905125BT	Pushrim Para 25" 6 bent tabs
art.no. 1905126BT	Pushrim Para 26" 6 bent tabs
art.no. 1905220BT	Pushrim Tetra 20" 6 bent tabs
art.no. 1905222BT	Pushrim Tetra 22" 6 bent tabs
art.no. 1905224BT	Pushrim Tetra 24" 6 bent tabs
art.no. 1905225BT	Pushrim Tetra 25" 6 bent tabs
art.no. 1905226BT	Pushrim Tetra 26" 6 bent tabs

NEW ARTICLE NUMBERS FOR THE NEW PUSHRIMS