

S3 CASTOR WHEEL ON PANTHERA X







Panthera S3 castor wheel



Panthera X with std castor

Information of what parts you need if you like to use the **new S3 castor wheel on Panthera X**.

These parts **MUST** be used when mounting the S3 castor wheel on Panthera X. Regardless if the wheelchair is on a new order or if you will change the castors on an existing wheelchair.



90 mm wheel			
Chassis	Rear wheel	art.no.	Description
X	24"	2100201	Fork 75, wheel 90 mm, axle X2/24"
X	25"	2100201	Fork 75, wheel 90 mm, axle X2/24"
X	26"	2100204	Fork 90, wheel 90 mm, axle X2/24"
120mm wheel			
Chassis	Rear wheel	art.no.	Description
X	26"	2100205	Fork 75, wheel 120 mm, axle X2/24"
120mm wheel on a 5" X chassis			
X 5"	24"	2100205	Fork 75, wheel 120 mm, axle X2/24"
X 5"	25"	2100205	Fork 75, wheel 120 mm, axle X2/24"
X 5"	26"	2100206	Fork 90, wheel 120 mm, axle X2/24"