


PREASSEMBLED CASTOR FORK FOR S3 AND X



S3 CASTOR FORK WITH AXLE



NOW WE DELIVER PART 1 AND 2 PREASSEMBLED

We will now deliver castor fork axle (1) and castor fork (2) preassembled from Panthera. We will not offer castor fork and axle separate. This we do to eliminate the risk of poorly mounted castor forks out on the field. This applies to S3 and X castor forks where the threaded axle is screwed to the castor fork and locked with thread locking fluid Loctite.

Preassembled S3 and X castor fork

Art.no: 2100140	Fork 75 +axle 3
Art.no: 2100141	Fork 90 +axle 3
Art.no: 2100142	Fork 105+axle 3
Art.no: 2100143	Fork 120+axle 3
Art.no. 2305139	Fork X2 axle 24" for X 3" X wheel
Art.no. 2305140	Fork X2 axle 25" for X 3" X wheel
Art.no. 2305141	Fork X2 axle 26" for X 3" X wheel