

## SEAT BAG IN A NEW DESIGN



SEAT MOUNTED art.no. 9900212



CUSHION MOUNTED art.no. 9900213

Now we release a seat bag in a new design. The seat bag is an accessory for storage of keys, cellphone or a wallet etc. It is delivered with the U3 LIGHT as standard.

The seat bag comes in two versions, seat mounted or cushion mounted.

-The **seat mounted** seat bag is attached directly to the velcro at the seat, and remain mounted to the wheelchair.

-The **cushion mounted** seat bag is attached to the cushion, and remain mounted to the cushion when lifted out of the wheelchair. The cushion mounted bag is delivered together with a self-adhesive velcro tape that you attached to the underside of the cushion. The velcro of the seat bag is then attached to this.

The seat bag's loose end is attached to the chassis tubes with magnets.

The seat bag is delivered together with U3 Light. In the orderguide for U3 Light you make a choice of what model of the seat bag you desire.

