

SPINERGY SPOX 24" EXPIRES AND IS REPLACED BY SPINERGY X 24"



SPINERGY SPOX 24"



SPINERGY X 24"

THE TRANSITION FROM SPINERGY SPOX TO SPINERGY X IS NOW COMPLETED WITH THE 24" WHEELS.

Soon we will start replacing the rear wheels Spinerger Spox 24" with Spinerger X 24", as we have done previously with Spinerger Spox 25" and 26".

Spinerger X is the standard rear wheel for Panthera X and it's a wheel with very low weight, high quality and stability.

The attachment of the push rim has the same design as all our other rear wheels which facilitate store-keeping. All our push rims fits all our rear wheels. Detachment and mounting of the pushrim is made from the outside, i.e. you don't need to remove the tire and tube for this operation.

The price for Spinerger X 24" will be the same as for Spinerger Spox 24"

Example of article numbers:

art.no. 1100081 Rear wheel 24" Spox std, without pushrim, Schwalbe Right Run tire

replaced by

art.no. 1100154 Rear wheel 24" Spinerger X, without pushrim, Schwalbe Right Run tire