

## PUSH HANDLES S3 H/S EXTENDED +10 CM

Clicka on image for larger format



PUSH HANDLES S3 H/S +10 CM

### WE NOW OFFER PUSH HANDLES S3 H/S EXTENDED +10 CM

Since we have introduced the upgraded backrest for the S3 series ( [see info sheet 2017-13](#)) we can now offer extended push handles S3 H/S +10 cm.

These push handles are extended 10 cm so that they can be adjusted to the right height for taller people, or when a lower backrest is mounted.

The extended push handles are only suitable for wheelchairs with the upgraded backrest VER2:

The push handles can be used with the following wheelchairs:

S3, S3 Short, S3 Short-Low, S3 large, S3 Swing, S3 Swing-Short och U3.

***(Note! Wheelchairs with a seat width of 45 cm and above have a maximum user weight of 150 kg. However, the push handle + 10 cm always has a maximum user weight of 100 kg)***

Art.no.	Description
4252063 VER2	Push handles S3 H/S +10 pair cpl
4252064	Push handles S3 H/S Left+10 (without attachment)
4252065	Push handles S3 H/S Right+10 (without attachment)