

## PRODUCT UPDATE

Click on images for larger format



CASTOR S2 5"



CASTOR S3 120 mm

### CASTOR WHEELS FOR THE S2 SERIES WILL STOP BEING PRODUCED

Castor wheels for the S2 series will stop being produced and will be replaced with castor wheel S3.

When mounting a S3 castor on a S2 castor fork you need "*Spacer S3-Casterwheel on SU2-fork*" **art.no. 9024002** (x2/fork)

We will offer these spacers at no cost when ordering spare parts until further notice. NOTE! You need to place an order on the spacers to get them delivered.

We will stop producing S2 castors and replace them with these S3 castors. We calculate that the change-over will be at the turn of the year 2018/2019 depending on quantity in stock.

The total cost for castors on a wheelchair calculated with a life span of 5 years will be, with this change, approx. 35 % lower with S3 castors compared to the S2 castors.

In this calculation we have considered the material cost and also the labor cost, that is; the new castors will have lower total cost despite a higher price for each castor.

EXPIRING		REPLACING	
Art.no.	Description	Art.no.	Description
2000040	Castor 3"/87mm STD	2002010	Castor S3 90 mm
2000000	Castor 5"/115mm standard rubber	2002020	Castor S3 120 mm
2000001	Castor 5"/120mm thick rubber	2002020	Castor S3 120 mm
2000032	Castor 6"/150mm	2002030	Castor S3 150 mm
2000041	Rubber thread 3"/87mm	2002010	Castor S3 90 mm
2000010	Rubber thread std 5"/115mm	2002020	Castor S3 120 mm
2000011	Rubber thread thick 5"/120mm	2002020	Castor S3 120 mm