

Product update

Click on images for larger format



WEDGE CUSHION



WEDGE CUSHION

Article: *Wedge cushion*

NEW PRODUCT: WEDGE CUSHION

Some customers have requested a flatter seat angle than our chairs currently have.

We can now offer a wedge cushion made from foam rubber designed for our wheelchairs.

The idea is for it to be placed under the usual seat cushion inside the cushion cover.

The wedge cushion is 3 cm high at the back, giving a flat seat angle for Swing and Bambino and halving the angle of the S3 and U3 models.

The depth is 41 cm and the cushion can be cut to the desired seat depth. [See image of measurement specs.](#)

DESCRIPTION	ART.NO.
Wedge cushion w. 27	5771027
Wedge cushion w. 30	5771030
Wedge cushion w. 33	5771033
Wedge cushion w. 36	5771036
Wedge cushion w. 39	5771039
Wedge cushion w. 42	5771042
Wedge cushion w. 45	5771045
Wedge cushion w. 50	5771050

For questions about this product contact:

Marketing manager

Nils Sallnäs

Tel: +46-8-564 73 862

E-mail: nils@panthera.se