

## U3 LIGHT BACKREST UPHOLSTERY FOR U2 LIGHT

Click on images for larger format



BACKREST UPHOLSTERY U3 LIGHT



BACKREST UPHOLSTERY U2 LIGHT

### BACKREST UPHOLSTERY FOR *U2 LIGHT* WILL BE REPLACED WITH *U3 LIGHT* UPHOLSTERY

We will replace the U2 Light upholstery with upholstery for U3 Light continuously when ordering spare parts, as they run out of stock.

The change will have logistic benefits and the U3 light upholstery has fewer seams against the body and better comfort. Slightly lower weight and a nicer design with a reflective band at the back.

Valid for art.no:

EXPIRES		REPLACES	
ART.NO.	DESCRIPTION	ART.NO.	DESCRIPTION
3263325LB	Upholstery backrest U2 light br 33 h 25	5263325L	Upholstery backrest kpl U3 light br 33 hō 25
3263330LB	Upholstery backrest U2 light br 33 h 30	5263330L	Upholstery backrest kpl U3 light br 33 hō 30
3263335LB	Upholstery backrest U2 light br 33 h 35	5263335L	Upholstery backrest kpl U3 light br 33 hō 35
3263625LB	Upholstery backrest U2 light br 36 h 25	5263625L	Upholstery backrest kpl U3 light br 36 hō 25
3263630LB	Upholstery backrest U2 light br 36 h 30	5263630L	Upholstery backrest kpl U3 light br 36 hō 30
3263635LB	Upholstery backrest U2 light br 36 h 35	5263635L	Upholstery backrest kpl U3 light br 36 hō 35
3263925LB	Upholstery backrest U2 light br 39 h 25	5263925L	Upholstery backrest kpl U3 light br 39 hō 25
3263930LB	Upholstery backrest U2 light br 39 h 30	5263930L	Upholstery backrest kpl U3 light br 39 hō 30
3263935LB	Upholstery backrest U2 light br 39 h 35	5263935L	Upholstery backrest kpl U3 light br 39 hō 35
3264225LB	Upholstery backrest U2 light br 42 h 25	5264225L	Upholstery backrest kpl U3 light br 42 hō 25
3264230LB	Upholstery backrest U2 light br 42 h 30	5264230L	Upholstery backrest kpl U3 light br 42 hō 30
3264235LB	Upholstery backrest U2 light br 42 h 35	5264235L	Upholstery backrest kpl U3 light br 42 hō 35
3264525LB	Upholstery backrest U2 light br 45 h 25	5264525L	Upholstery backrest kpl U3 light br 45 hō 25
3264530LB	Upholstery backrest U2 light br 45 h 30	5264530L	Upholstery backrest kpl U3 light br 45 hō 30
3264535LB	Upholstery backrest U2 light br 45 h 35	5264535L	Upholstery backrest kpl U3 light br 45 hō 35