

NEW SOLID PUNCTURE FREE PU- TIRE

Click on images for larger format



PREVIOUS PU -TIRE

NEW PR1MO PU -TIRE

NEW SOLID PUNCTURE FREE PU- TIRE

Starting in August 2019, we will continuously replace the old gray PU tires (solid polyurethane puncture free) with our new black PR1MO PU tires.

The change means that we will have a solid tire with better braking action than previous version, it's quicker to mount in a tire mounting - jig, due to it's easier to fit it into the rim. In addition, the new tire rolls substantially better, according to our tests about 30% less rolling resistance. Will also be available in size 25 "and 26".

Available in sizes 20", 22", 24", 25" and 26".

The PR1MO tire fits both our [std wheels](#) and the [Spinergywheels](#).

NOTE! Does not fit the rear wheel for [one arm drive](#) or the [assistant brake](#). For those we use the [gray PU tire](#).

The price remains the same. Article number for the std rearwheel with premounted PU-tire will also be unchanged.

Assembly video: [VIDEO](#)

Valid for art.nr:

EXPIRES		REPLACES			REARWHEELS WITH PU-TIRE
ART.NO.	DESCRIPTION	ART.NO.	DESCRIPTION	ART.NO.	DESCRIPTION
1000116	Tire solid 20"	1000220	Tire solid 20"x1" Black Primo	1320001	Rearwheel, 20"6P no pushrim solid PU tire
1000117	Tire solid 22"	1000222	Tire solid 22"x1" Black Primo	1310001	Rearwheel, 22"6P no pushrim solid PU tire
1000115	Tire solid 24"x1"	1000224	Tire solid 24"x1" Black Primo	1300001	Rearwheel, 24"6P no pushrim solid PU tire
		NEW			
		1000225	Tire solid 25"x1" Black Primo		
		1000226	Tire solid 26"x1" Black Primo		