

## New bushing for the S3 castor wheel

Click on images for larger format



S3 CASTOR WHEEL ON FORK

CASTOR WHEEL BUSHINGS

PREVIOUS CASTOR WHEEL BUSHINGS

NEW CASTOR WHEEL BUSHING

UPCOMING CASTOR FORK S3

MOUNTING WITH NEW FORK

We now introduce a new design of the bushing for the S3 castor wheel. This is to simplify the logistics and the mounting of the castor wheel onto the castor fork.

The previous bushings for the S3 castor wheels came in two versions, one for the left and one for the right side. Now we replace these two with one version that fits both sides. With only one edge on the new bushing it becomes easier to mount.

We will also later on, introduce an updated castor fork. When the lower edge disappears, [see image](#), it becomes much easier to mount the castor wheel onto the castor fork.

Both versions of the bushings fits both versions of the castor fork.

ART.NO.	DESCRIPTION	
2002017	Bushing casterwheel S3	New version, for both left and right side
2002014	Bushing casterwheel S3, right	Previous version, for the right side
2002015	Bushing casterwheel S3, left	Previous version, for the left side