

## QR axle for rear wheel Spinerger X 24" and 25"

Click on image for larger format



QR axle for rear wheel Spinerger X

X, left - U3 light, right

Together with rear wheel Spinerger X in 24" and 25", we deliver the QR axle 1/2" 12,7 mm x 99 mm titanium with art.no. 1700710.

During an introduction period for the wheelchair model U3 Light, this axle was delivered with a U3L label.

This label will now expire and all rear wheel axles delivered will have the original marking with the "X" logotype.

ART.NO.	DESCRIPTION
1700710	QR axle 1/2" 12,7mm x 99mm titanium