

Product improvement S3 Anti-tip device for low wheelchairs

Click on image for larger format



ANTI-TIP DEVICE S3

S3 Anti-tip device for low wheelchairs

Until now, you have had to modify the anti tip devices S3 to be mounted on our low wheelchair models S3 Short Low, S3 Swing Low 0°, U3 Y-front Low.

The reason is that the position of the lock pin on these models and it has been necessary to change the left and right sides during assembly and drill holes for the knob.

Now we can offer Anti tip device S3 Low where you don't have to adapt before installation.

The price is the same as for [Anti tip device S3](#)

New article numbers :

- 4355000 Anti tip device S3 Low cpl pair
- 4355001 Anti tip device S3 Low cpl left
- 4355002 Anti tip device S3 Low cpl right
- 4353430 Anti tip device S3 Low +30 pair
- 4353431 Anti tip device S3 Low +30 left
- 4353432 Anti tip device S3 Low +30 right

Please do not hesitate with any questions:

Export / Product Manager

Helena Lemoine

helena.lemoine@panthera.se

Tel: +46-8-557 608 77