

To: Prescribers
Technicians
Purchasers

Backward bent back frame on S2 swing

Due to the 4 degrees seat-angle on the S2 Swing the special backrest that's curved/bent backward is not possible to angle backwards as much as on all the other frame-models with 7 degree seat-angle. (Its OK for approx. 60% of the users)

We are changing the design of the backrest but we will not have any alternative until earliest June.

If you anyway order this combination and the customer is not satisfied you can replace the backrest-tube as soon as we have a substitute.

Backward curved
frame

