

Adjusted by:

E-mail:

Adjusted for:

Tel:

Date:

Total weight (g):

The weight changes based on the choices made below.

S3 short has a adjustable rear axle in carbon fibre. To find a balance that suits you the rear axle can be adjusted into different positions. Note desired measure at "Settings". S3 short has a seat height of 47 cm in the front and 43 at the rear, without cushion. To adjust the seat height you can chose other rearwheel sizes or a raising Kit with suitable castor fork from the List of options.

Max userweight: 100 kg for seat width 30-42 and 150kg for seat width 45 cm

Save the guide to the "desktop" and open it with [Adobe Reader](#), for best function Design your wheelchair by clicking in the circles in this guide. Blue, underlined text is clickable to [show info](#) or image. BB = seatwidth, HH= backrest height

Original design

Alternative options

Other options

[info](#)

1) Chassis-colour

- Black - Anodic Black
- Gray-Gray metallic
- Clear

- [Optional colour](#) - A higher price and longer delivery time for optional colour, Fili in art.no. in the square



2) Seat width (BB)

- [Art.no.](#)
- Clear
- 30 cm
- 33 cm
- 36 cm
- 39 cm
- 42 cm
- 45 cm

Art.no include: chassis, seat, calf strap, backrest lock

3) Seat depth

- Clear
- Adjustable seat depth 27-33 cm
Standard for seat width 30-33 cm
- Seat depth 35 cm
For all seat width



4) Backrest

- [Standard](#)
523BBHH VER2
 - Clear
- Height (HH)
- 20 cm
 - 25 cm
 - 30 cm
 - 35 cm
 - 40 cm
 - 45 cm
 - Clear

- [Tapered](#)
525BBHH VER2
- [Curved back](#)
524BBHH VER2

[Other backrest](#)

[List of options](#) art. _____



5) Rear axle

- [Carbon standard](#)
33023BB
- Clear

- [Carbon +15 mm](#)
33023BC*
- * seat width BB + 1
- Use this rear axle when chosing arm support at Accessories and/or increased seat height

[Other](#) rear axle or seat height

[Increase the seat height](#)

[List of options](#) art. _____



6) Castor and fork

- [120 mm](#)
[art.no.](#)
- Clear

- [90 mm](#)
[art.no.](#)
- [150 mm](#)
[art.no.](#)



7) Footrest

- [Foorest Titanium](#)
36001BB
- [Footrest U-shape small heel bar](#)
3600131 (width 30) 3600134 (width 33)
- Clear

- [Extended, reinforced](#)
36111 BB [needed for Freewheel](#)
- [Foldable footplates](#)
[Art.no.](#)
- [Heelstrap for footplates](#)
[Art.no.](#)
- [Plastic plate](#)
36010BB

[Other footrest](#)

[List of options](#) art. _____



8) Brake

- [Brake S3](#)
4630000
- Clear

- [One hand brake](#)
[Art.no.](#)

[Other brake](#)

[List of options](#) art. _____



9) Rear wheel - without pushrim

- 24" [standard](#)
[art.no.](#)
- Clear

- 22" [standard](#) [Rearwheel sizes](#)
[art.no.](#)
- 25" [Spinergy X](#) black spokes
[art.no.](#)
- 26" [Spinergy X](#) black spokes
[art.no.](#)

[Other rear wheel](#)

[List of options](#) art. _____



Tire - choose tire

- [Schwalbe Right Run](#) 22",24",25",26"
- Clear

- [Rugged](#) 24"
- Schwalbe [Marathon Plus](#) 24",25",26"

[PU](#) 24" (massive)

[reTyre](#)-tyre skins (does not fit PU tyres)



10) Pushrim

- [Titanium](#)
190DD11 (x2 / wheelchair)
 - Clear
- DD = Wheel diameter

- [Aluminum](#)
1904404 x2 (04=24", 05=25", 06=26")
- [Paragrip](#)
19051DDBT x2
- [Tetragrip](#)
19052DDBT x2

[Maxgrip](#) 22", 24" (contains latex)
1903060 x2 (60=24", 80=22")

Accessories



Side guard

- [Side guard S3](#) [Video](#)
[art.no.](#)
- Clear



Arm support

- [Arm support S3](#) [Video](#)
[art.no.](#)



Side guard Carbon

- [Side guard Carbon](#)
[art.no.](#)



Push handles

- [Height adjustable S3](#) [Video](#)
4252060 VER2
- [Height adjustable S3 +10](#)
4252063 VER2
Max. userweight 100kg
- Clear



Push handle Foldable

- [Foldable S3](#) [Video](#)
4245000



Push bar

- [Push bar S3](#)
42400BB VER2
Max. userweight 100kg



Anti-tip


- [Anti-tip](#) std seat height
4353000 [Video](#)

[Other seat height](#)



Cushion Wedge cushion

- [2,5 cm](#)
572BB35
- [5 cm](#)
575BB35
- Clear
- [Wedge cushion](#)
57710BB



Spoke protectors

- [24" logo](#)
4500000
- [22" logo](#)
4500100
- [24" no logo](#)
4500011
- [22" no logo](#)
4500111
- Clear



Other

- [Tube protection front](#)
9900251 (2 pcs)
- [Tube protection rear](#)
9900251 (1 pcs)
- [Crutch holder S3](#) _____ pcs
- [9900301](#) (x1/crutch)
- [Backrest wedge](#) _____ pcs
4710012

- [Table S3](#)
6810020



Chose art.no. for Freewheel according to the **footrest height** of the wheelchair
Read more about Freewheel on the [homepage](#)

- [Freewheel](#)

Set-up adjustments



Back rest angle

Note measure (A) for desired back rest angle.
14-28 mm. Standard=20 mm

A= _____ mm



Rear axle adjustment

Note measure (A) for desired balance.
Standard=40 mm


A= _____ mm



Footrest height **Footrest - Foot plates**

Note measure (A) and (B) for desired footrest.
[Bambino foot plate with Feewheel - info](#)

A= _____ mm B= _____ mm



Transportation marking

YES do mark the wheelchair with
 attachment points for transport in
disabled service vehicle. Art.no.9000152



NOTE! See [manual](#) for more information.

Extra