

Adjusted by:

E-mail:

Adjusted for:

Tel:

Date:

Total weight (g):

Vikten ändrar sig beroende på vilka val du gör nedan.

U3-light Y-front 90° has a rigid rear axle. To find a balance that suits you, there are two versions, **U3 light Y-front 90°** or **U3 light Y-front L 90°**. L is more easy balanced and easier to drive on the rear wheels. You can tune the balance by placing the backrest in three different position. U3 Light has a seat height of 47 cm in front and 43 cm at the rear without cushion. **Max userweight is 100 kg for seatwidth 33-45.**

Save the guide to the "desktop" and open it with **Adobe Reader**, for best function **Design your wheelchair by clicking in the circles in this guide.** **Blue, underlined text is clickable to show info** or image. **BB** = seatwidth, **HH** = backrest height

Original design

Alternative options

Other options

[info](#)

1) Chassis-colour Click on the link to see up-to-date colours.

Fill in art.no. in the square



2) Seat width (BB)

[art.no.](#)

33 cm

36 cm

39 cm

42 cm

45 cm

Clear

Art.no include: chassis, seat, calf strap, and seat bag



3) Seat depth

[art.no.](#)

Rensa

35 cm

37,5 cm

40 cm

[info 40 cm](#)



4) Back rest

Height (HH)

U3 light

523BBHH L

Clear

25 cm

30 cm

35 cm

Clear

Tapered

525BBHH L

Bakåttbockad

524BBHH L

Upholstery, comfort

U3 Light backrest with S3 upholstery - more padding and stronger ribbons than std - slightly heavier

Backrest, extended range

Frame fo fixed backrest systems

art. _____



5) Balance center of gravity (COG)

[info](#)

Position A
Least easy balanced

Position B
Medium balance

Position C
Most easily balanced

Rensa

The positioning of the backrest on the chassis affects the balance / COG of the wheelchair. "easy balanced" = easy to practice rear wheel balance. This can later easily be changed.



6) Castor and fork

90 mm X

[art.no.](#)

Rensa

90 mm S3

[art.no.](#)

120 mm S3

[art.no.](#)



7) Footrest

Footrest Y3 90

36018BB

Clear

Plastic plate for footrest

36010BB



8) Brake

One arm brake, right

[art.no.](#)

Clear

Brake S3

4630000

One arm brake, left

[art.no.](#)



9) Rear wheel - without pushrim [art.no.](#)

24" Spinerger X white spokes

24" Spinerger X black spokes

Clear

25" Spinerger X white spokes [Drivhjulsstorlek](#)

25" Spinerger X black spokes

26" Spinerger X white spokes

26" Spinerger X black spokes

Other rear wheel



Tire - choose tire

Schwalbe One 24", 25"

Clear

Right Run 24", 25" 26"

Marathon Plus 24", 25", 26"

Rugged 24"

reTyre-tyre skins



10) Pushrim

- [Aluminum X 24"](#)
1904404 (x2)
- [Aluminum X 25"](#)
1904405 (x2)
- [Aluminum X 26"](#)
1904406 (x2)
- Clear

- [Titanium](#)
190DD11 x2
- [Paragrip](#)
19051DDBT x2
- [Tetragrip](#)
19052DDBT x2

DD = Hjulldiameter

- [Maxgrip 24"](#) (contains latex)
1903060 x2

Accessories



Side guard Carbon

- [Carbon](#)
art.no.
- Clear



Side guard

- [Carbon, small cover](#)
art.no.
 - [Carbon, large cover*](#)
art.no.
 - [Side guard S3](#)
art.no.
- * Only with one arm brake



Arm support

- [Arm support S3](#) [Video](#)
art.no.

Makes the wheelchair 15mm wider



Cushion


Seat depth 35cm Seat depth 40cm

- [2,5 cm](#) [2,5 cm](#)
572BB35 572BB40
- [5 cm](#) [5 cm](#)
575BB35 575BB40
- Clear



Push handle Foldable


- [Foldable S3](#) [Video](#)
4245000
- Clear



Seat bag

Included in the prize. Choose model.

- [Seat mounted](#)
9900212
- [Cushion mounted](#)
9900213
- Clear



Wedge cushion


- [Wedge cushion](#)
57710BB



Other

- [Tube protection](#) ___ p
9900251 (x3/wheelchair)
- [Backrest wedge](#) ___ p
4710012 Chose art.no. for Freewheel according to the **footrest height** of the wheelchair Read more about Freewheel on the [homepage](#)
- [Freewheel](#)

Set-up adjustments



Back rest angle

Note measure (A) for desired back rest angle.
10-20 mm. Standard=18 mm

A= _____ mm



Footrest height [info about footrest height](#)

Ange måttet (A) för önskad fotstödshöjd.

A= _____ mm

Extra