

Adjusted by:

E-mail:

Adjusted for:

Tel:

Date:

Total weight (g):

The weight changes based on the choices made below.

U3-light Y-front Low has a rigid rear axle. To find a balance that suits you, there are two versions, **U3 light Y-front or U3 light Y-front Low L**. **L** is more easy balanced and easier to drive on the rear wheels. You can tune the balance by placing the backrest in three different position. U3 Light has a seat height of 44,5 cm in front and 41,5 cm at the rear without cushion.
Max userweight is 100 kg for seatwidth 33-45.

Save the guide to the "desktop" and open it with **Adobe Reader**, for best function Design your wheelchair by clicking in the circles in this guide. **Blue, underlined text is clickable to show info** or image. **BB** = seatwidth, **HH** = backrest height

Original design

Alternative options

Other options

[info](#)

1) Chassis-colour

Click on the link to see up-to-date colours.

Fill in art.no. in the square



2) Seat width (BB)

[art.no.](#)

- 33 cm
- 33 cm L
- Clear

- 36 cm
- 36 cm L

- 39 cm
- 39 cm L

- 42 cm
- 42 cm L

- 45 cm
- 45 cm L

Art.no include: chassis, seat, calf strap, and seat bag



3) Seat depth

[art.no.](#)

- Clear
- 35 cm
- 40 cm

- 37,5 cm
- 42,5 cm



4) Back rest

Height (HH)

- U3 light
523BBHH L
- Rensa
- 25 cm
- 30 cm
- 35 cm
- Rensa

- Tapered
525BBHH L
- Curved backd
524BBHH L

[Info](#)

- Upholstery, comfort
U3 Light backrest with S3 upholstery - more padding and stronger ribbons than std - slightly heavier

- Backrest, extended range
- Frame fo fixed backrest systems

art. _____



5) Balance center of gravity (COG)

[info](#)

- Position A
Least easy balanced
- Position B
Medium balance
- Position C
Most easily balanced
- Clear

The positioning of the backrest on the chassis affects the balance / COG of the wheelchair. "easy balanced" = easy to practice rear wheel balance. This can later easily be changed.



6) Castor and fork

- 90 mm X
[art.no.](#)
- Clear

- 90 mm S3
[art.no.](#)
- 120 mm S3
[art.no.](#)



7) Footrest [info about footrest height](#)

- Footrest X
36004BB
- Clear

- Footrest Titanium
36915BB
- Extended, reinforced
36006BB [needed for Freewheel](#)
[Plastic plate for footrest](#)
36010BB

- Other footrest



8) Brake

- One arm brake, right
[art.no.](#)
- Clear

- Brake S3
4630000

- One arm brake, left
[art.no.](#)



9) Rear wheel - without pushrim [art.no.](#)

- 24" Spinerger X white spokes
- 24" Spinerger X black spokes
- Clear

- 25" Spinerger X white spokes [Drivhjulstorlek](#)
- 25" Spinerger X black spokes
- 26" Spinerger X white spokes
- 26" Spinerger X black spokes

- Other rear wheel



Tire - choose tire

- [Schwalbe One 24", 25"](#)
- Clear

- Right Run 24", 25" 26"
- Marathon Plus 24", 25", 26"

- Rugged 24"

- reTyre-tyre skins (does not fit PU tyres)



10) Pushrim

- [Aluminum X 24"](#)
1904404 (x2)
- [Aluminum X 25"](#)
1904405 (x2)
- [Aluminum X 26"](#)
1904406 (x2)
- Clear

- [Titanium](#)
190DD11 x2
- [Paragrip](#)
19051DDBT x2
- [Tetragrip](#)
19052DDBT x2

DD = Wheel diameter

- [Maxgrip 24"](#) (contains latex)
1903060 x2

Accessories



Side guard *

- [Carbon, without cover](#)
[art.no.](#)
- * NOTE! Side guards for 26" is required for 24" rear wheels. Due to the low chassis
- Clear



Side guard *

- [Carbon,small cover](#)
[art.no.](#)
- [Carbon, large cover*](#)
[art.no.](#)
- * Only with one arm brake
- [Side guard S3](#)
[art.no.](#)



Arm support *

- [Arm support S3](#) [Video](#)
[art.no.](#)
- * Makes the wheelchair 15mm wider
Not possible to combine with 25" and 26" rear wheels.
Arm support 26" is required for 24" rear wheels due to the low chassis



Cushion

Seat depth 35cm Seat depth 40cm

- [2,5 cm](#) [2,5 cm](#)
572BB35 572BB40
- [5 cm](#) [5 cm](#)
575BB35 575BB40
- Clear



Push handle Foldable


- [Foldable S3](#) [Video](#)
4245000
- Clear



Seat bag

Included in the prize. Choose model.

- [Seat mounted](#)
9900212
- [Cushion mounted](#)
9900213
- Clear



Wedge cushion


- [Wedge cushion](#)
57710BB



Other

- [Tube protection](#) ___ p
9900251 (x3/wheelchair)
- [Backrest wedge](#) ___ p
4710012
- Chose art.no. for Freewheel according to the **footrest height** of the wheelchair Read more about Freewheel on the [homepage](#)
- [Freewheel](#)

Set-up adjustments



Back rest angle

Note measure (A) for desired back rest angle.
10-20 mm. Standard=18 mm

A= _____ mm



Footrest height [info about footrest height](#)

Note measure (A) for desired footrest height.

A= _____ mm

Extra