

Adjusted by:
 E-mail:
 Adjusted for:

Tel:
 Date:

Total weight (lb):

The weight changes based on the choices made below.

U3 has a adjustable rear axle in carbon fibre. To find a balance that suits you the rear axle can be adjusted into different positions. Note desired measure at "Settings". U3 has a seat height of 18,5" in the front and 16,9" at the rear, without cushion. To adjust the seat height you can chose other rearwheel sizes or a raising Kit with suitable castor fork from the List of options.
Max userweight: 220 lb for seat width 13-13,5" and 330 lb for seat width 18"

Save the guide to the "desktop" and open it with **Adobe Reader**, for best function
Design your wheelchair by clicking in the circles in this guide.
 Blue, underlined text is clickable to [show info](#) or image.
BB = seatwidth(cm), **HH**= backrest height (cm) [INFO about measures](#)

Original design

Alternative options

Other options

[info](#)

1) Chassis-colour

- Black** - Anodic Black

[INFO about measures](#)

- Optional colour - A higher price and longer delivery time for optional colour

[Optional colour](#)

2) Seat width (BB in cm)

[article number](#)

- Clear

- 13"
- 14"
- 15"
- 16,5"
- 18"

Art.no include:
chassis, seat, calf strap, backrest lock

3) Back rest

- Standard**
523BBHH

- Clear

Height (HH)

- 8"
- 10"
- 12"
- 14"
- 16"
- 18"

Clear

- Tapered**

525BBHH

[Info](#)

- Curved back**

524BBHH

- Other back rest**

[List of options](#) art. _____

4) Rear axle

- Carbon standard**
33023BB

- Clear

- Carbon +0,59"**

33023BC*

* seat width BB +1

Use this rear axle when choosing arm support at Accessories

- Other** rear axle or seat height

[Increase the seat height](#)

[List of options](#) art. _____

5) Castor and fork

- 3" /90 mm**
2100030 x2

- Clear

- 5" / 120 mm**

2100020 x2

6) Footrest

- Footrest Titanium**
36915BB

- Clear

- Extended, reinforced**

36006BB

- Plastic plate for footrest**

36010BB

- Other footrest**

[List of options](#) art. _____

7) Brake

- Brake S3**
4630000

- Clear

- One hand brake**

[Art.no.](#)

- Other brake**

[List of options](#) art. _____

8) Rear wheel - without pushrim

- 24" **standard** with std tire
1300000 (x2 / wheelchair)

- Clear

- 24" **Spinergy Spox**

1100081 x2

- 25" **Spinergy Spox**

1100130 x2

- 26" **Spinergy Spox**

1100100 x2

[Rear wheel sizes](#)

- Other rearwheels**

[List of options](#) art. _____

9) Push rim

- 24" **Titanium**
190DD11 (x2 / wheelchair)

- Clear

DD = Wheeldiamter

- Aluminum**

190DD30 x2

- Tetragrip**

19052DDBTx2

- Paragrip**

19051DDBTx2

- Maxgrip** (contains latex)
1903060 x2

10) Tire

- Schwalbe Right Run** 24",25",26"
Standard


- Clear

[art.no.](#)

- Rugged** 24"

- Schwalbe **Marathon Plus**
24",25",26"

- PU** 24" (massive)



Accessories



Side guards

- [Side guards S3](#) [Video](#)
article number
- Clear



Arm support

- [Arm support S3](#) [Video](#)
article number



Sidoskydd U2 light

- [Side guards U2 light](#)
article number



Push handles

- [Height adjustable S3](#) [Video](#)
4252060
- Clear



Push handle Foldable

- [Foldable S3](#) [Video](#)
4245000



Push handle Pushbar

- [Pushbar S3](#)
42400BB



Anti-tip

- [Anti-tip](#) std seat height
4353000 [Video](#)

[Other seat height](#)



Cushion

- [1" / 2.5 cm](#)
47200BB
- [2" / 5 cm](#)
47500BB
- Clear



Spoke protectors

- [24" logo](#) [22" logo](#)
- 4500000 4500100
- [24" no logo](#) [22" no logo](#)
- 4500011 4500111
- Clear



Other

- [Tube protection](#) _____ pcs
- 9900251 (x3/wheelchair)
- [Crutch holder](#) _____ pcs
- 9900300 (x2/crutch)
- [Back rest wedge](#) _____ pcs
- 4710010

Set-up adjustments



Back rest angle

Note measure (A) for desired back rest angle.
0.55"-1.1". Standard=0.8"

A= _____ "



Rear axle adjustment

Note measure (A) for desired balance.
Standard=2,36"

A= _____ "



Footrest height

Note measure (A) for the footrest.

A= _____ "



Transportation marking

- YES do mark the wheelchair with attachment points for transport in disabled service vehicle

NOTE! See the [manual](#) for more info.

Extra