

Adjusted by:

E-mail:

Adjusted for:

Tel:

Date:

Total weight (lb):

The weight changes based on the choices made below.

**U3-light** has a rigid rear axle. To find a balance that suits you, there are two versions, **U3 light** or **U3 light L**. L is more easy balanced and easier to drive on the rear wheels. You can tune the balance by placing the backrest in three different position. U3 Light seat height of 18,5" in the front and 16,9" at the rear without cushion. **Max userweight is 220 lb for seatwidth 13-18"**.

Save the guide to the "desktop" and open it with **Adobe Reader**, for best function Design your wheelchair by clicking in the circles in this guide. **Blue, underlined text is clickable to show info** or image. **BB** = seatwidth (cm), **HH**= backrest height (cm) **INFO about measures**

**Original design**

**Alternative options**

**Other options**

[info](#)

**1) Chassis colour**

- White** - Pearl White  Clear
- Black** - Anodic Black

[INFO about measures](#)

**Optional colour** - A higher price and longer delivery time for optional colour

[Optional colour](#)



**2) Seat width (BB)**

[article number](#)

- 13"  14"  15"  16,5"  18"
- 13" L  14" L  15" L  16,5" L  18" L
- Clear

Art.no includes: chassis, seat, calf strap and a seat bag



**3) Backrest**

Height (HH)

- U3 light**  10"
- 523BBHH L  12"
- Clear  14"
- Clear

**Backrest upholstery Comfort**

- U3 Light backrest with S3 upholstery - more padding and stronger ribbons than std - slightly heavier



**4) Balance center of gravity (COG)**

[info](#)

- Position A**  **Position B**  **Position C**
- Least easy balanced Medium balance Most easily balanced
- Clear

The positioning of the backrest on the chassis affects the balance / COG of the wheelchair. "easy balanced" = easy to practice rear wheel balance. This can later easily be changed.



**5) Castor**

- 3" / 90 mm X**
- 2100222 (x2 / wheelchair)
- Clear

- 3"/90 mm S3**
- 2100030 x2
- 5" / 120 mm S3**
- 2100020 x2



**6) Footrest**

- Footrest X**
- 36004BB
- Clear

- Footrest, Titanium**
- 36915BB
- Footrest extended, reinforced**
- 3606BB **needed for Freewheel**
- Plastic plate for foot bow**
- 36010BB

**Other footrest**

[List of options art.](#) \_\_\_\_\_



**7) Brake**

- One hand brake, right**
- [article number](#)
- One hand brake, left**
- Clear

- Brake S3**
- 4630000



**8) Rear wheel - without pushrim**

- 24" Spinergy X white spokes**
- 1100184 (x2 / wheelchair)
- 24" Spinergy X black spokes**
- 1100141 (x2 / wheelchair)
- Clear

- 25" Spinergy X white spokes**
  - 1100185 x2
  - 25" Spinergy X black spokes**
  - 1100146 x2
  - 26" Spinergy Spox black spokes**
  - 1100100 x2
- [Rear wheel sizes](#)

**Other rear wheel**

[List of options art.](#) \_\_\_\_\_



**9) Pushrim**

- Aluminum X**
- 1904404 (5) (x2 / rullstol) 24" (25")
- Clear

- Titanium**
- 190DD11 x2
- Paragrip**
- 19051DDBTx2
- Tetragrip** **DD = Wheel diameter**
- 19052DDBTx2

**Maxgrip 24"** (contains latex)

1903060 x2



**10) Tire**

- Schwalbe One 24", 25"**
- standard [art.no.](#)
- Clear

- Schwalbe Right Run**
- 24",25",26"
- Rugged 24"**

Schwalbe **Marathon Plus**

24",25",26"

## Accessories



### Side guard

- [U2 light](#) [info](#)  
article number
- [Side guard S3](#)  
article number
- Clear



### Side guard

- [Carbon, small cover](#)  
article number
- [Carbon, large cover](#)  
article number



### Arm support

- [Arm support S3](#) [Video](#)  
article number

Makes the wheelchair 15mm wider



### Cushion

- [1" / 2,5 cm](#)  
572BB40
- [2" / 5 cm](#)  
575BB40
- Clear



### Push handles

- [Foldable S3](#) [Video](#)  
4245000
- Clear



### Other

- [Freewheel info](#)  
2600102
- [Tube protection](#)  
9900251 (x3 / wheelchair) \_\_\_ pcs
- [Back rest wedge](#)  
4710010 \_\_\_ pcs

## Set-up adjustments



### Backrest angle

Note measure (A) for desired back rest angle.  
0,39-0,79 " Standard=0,71"

A= \_\_\_\_\_ "



### Footrest height

Note measure (A) for desired footrest height.  
Standard = 2,75"

A= \_\_\_\_\_ "

Extra